



Member Tested and Recommended Program

Tuesday, June 13, 2006

Thank you for participating in our “Member Tested and Recommended Products Program.” The **Opti Fit CBRN Mask** has been tested and recommended by the members of the National Tactical Officers Association. Your overall score was a **4.77**. The NTOA Member Tested and Recommended logo may be used solely on the specific product that was tested.

Attached with this letter is the Member Tested Logo, which can be used on marketing pieces for this product and a copy of the reviews. This review will soon be published in the Tactical Edge journal, the Member Tested Product Review Guide, and NTOA website.

If your company is interested in participating in this comprehensive review guide, please contact us. Once again thank you for participating in this worthwhile program.

Thank you,
National Tactical Officers Association



NTOA PO Box 797, Doylestown, PA 18901 www.ntoa.org

Defense Technology® Opti Fit CBRN Mask — Overall Score: 4.77

<p>TESTED BY A POLICE OFFICER FROM KENTUCKY</p> <p>Design 5 Performance 5 Ease of Use 5 Size 5 Quality 5 Durability 5 Storage 5 Versatility Convenience Application 5 Comfort 5 Accuracy Cleaning & Maintenance 5 INDIVIDUAL SCORE 5</p>	<p>I tested the Defense Technology Opti Fit CBRN (Chemical, Biological, Radiological and Nuclear) mask from the standpoint of comfort and application. The Opti Fit mask is really just as comfortable as any other mask that is on the market. What I really like about this mask is it offers the operator a greater field of view along the sides and bottom. The drinking tube worked well with the standard Camelback gear and all filters, tubes, etc. seemed to work just as well.</p> <p>I didn't get a chance to test it in an environment like OC. However, I used the standard test of placing my palm over the filter and inhaled. I also did some jumping-jacks, push-ups, and a little running with the mask, helmet and vest. The mask stayed in place and did not interfere with any of my other gear.</p>
---	--

<p>TESTED BY A POLICE OFFICER FROM NEW MEXICO</p> <p>Design 5.0 Performance 5.0 Ease of Use 4.0 Size 4.5 Quality 5.0 Durability 4.0 Storage 4.0 Versatility 5.0 Convenience 4.0 Application 5.0 Comfort 5.0 Accuracy Cleaning & Maintenance 4.0 INDIVIDUAL SCORE 4.54</p>	<p>The SURVIVAIR Opti-Fit Tactical Gas Mask was tested in a gas chamber and at the shooting range. The fit of the SURVIVAIR Opti-Fit Tactical Gas Mask was very comfortable and offered great unobstructed vision and no distortion. The Opti-Fit was tested with shoulder fired weapons utilizing different optical devices, (Eotech, Aim Point and CQT-4). No observable distortion was noticed while using the different optics with the SURVIVAIR Opti-Fit Tactical Gas Mask and breathing resistance was minimal. The construction of the full face-piece is soft silicone and was very comfortable during extended wear with a standard nose-cup to help reduce fogging. The SURVIVAIR Opti-Fit Tactical Gas Mask offers three positions to attach a 1688 NBC Canister, left side or right side for individual shooters and a conventional front mounting position. Third party testing of the SURVIVAIR Opti-Fit Tactical Gas Mask by the US Army Soldier Biological Chemical Command (SBCCOM) against HCN, CN, CS, Phosgene, Sarin and Mustard provides assurance that this respirator is effective against chemical warfare agents. I would recommend the SURVIVAIR Opti-Fit Tactical Gas Mask for any law enforcement agency/officer that may be looking for a great quality APR to operate in a CBRN environment.</p>
--	--